

24-2 & 24-2C Visual Field Exams

Your Comprehensive Guide

[24-2 Visual Field Exam Overview](#)

The 24-2 Visual Field Exam at a Glance

List of Ocular Diseases Monitored and Diagnoses Identified by the 24-2 Exam

Example 24-2 Reports

The 24-2C Visual Field Exam History

Billing and Coding for the 24-2 Eye Exam

Start Conducting the 24-2 Exam with Carrot

Screening patients for glaucoma or neuro-ophthalmic conditions often starts with the 24-2 visual field exam. This cornerstone diagnostic test systematically evaluates the central 24 degrees of a patient's vision and reveals patterns of visual field loss that may otherwise go unnoticed. You can use the test to achieve an initial diagnosis or track disease progression for glaucoma, optic neuropathies, brain lesions, and many other conditions that affect the visual pathways.

This guide simplifies the essentials: how the test works, its clinical applications, and how to use Carrot to upgrade the experience to improve patient comfort, precision, and ease.

24-2 Visual Field Exam Overview

In ophthalmology eye exams, the 24-2 visual field test assesses a patient's peripheral vision and detects abnormalities. It divides the central 24 degrees of the field of vision into a grid of 54 points. Each point is individually tested for light sensitivity as lights adjust in intensity to determine the patient's threshold at each location. Test results map the patient's visual field profile into graphs and patterns you can interpret for diagnosis.

It's best-known as a glaucoma detection test that reliably identifies patterns like nasal steps and arcuate scotomas, but it can diagnose and monitor many more conditions.

Occasionally, you'll see the 24-2 test called the HVF 24-2 or the [Humphrey visual field exam](#). This name is a reference to the Humphrey visual field perimeter, which has been historically used to conduct this test since the 1980s. However, using a Humphrey Visual Field Analyzer to conduct this exam is not strictly necessary. Carrot is an innovative, streamlined, and patient-friendly alternative.

[24-2 Visual Field Exam Overview](#)

The 24-2 Visual Field Exam at a Glance

List of Ocular Diseases Monitored and Diagnoses Identified by the 24-2 Exam

Example 24-2 Reports

The 24-2C Visual Field Exam History

Billing and Coding for the 24-2 Eye Exam

Start Conducting the 24-2 Exam with Carrot

Academic references and clinical validation



[A prospective cohort study](#) demonstrated that baseline central visual field loss measured with the 24-2 visual field testing strategy was predictive of more rapid and significant progression of global field damage in eyes with manifest glaucoma. This finding supports the clinical relevance and predictive value of the 24-2 exam in glaucoma management.



[Research by Hood et al. \(2022\)](#) compared the diagnostic accuracy of 24-2 and 10-2 perimetry tests. They found that the 24-2 test had better overall diagnostic accuracy, with a higher area under the ROC curve (0.808) compared to the 10-2 test (0.742). Results support the validity of the 24-2 exam by demonstrating its superior diagnostic accuracy compared to another commonly used test, reinforcing its effectiveness in clinical practice.



[Recent research by Kahook and Noecker](#) highlighted the importance of proper interpretation of the 24-2 Humphrey Visual Field printout, emphasizing factors such as pupil diameter and appropriate refraction that can influence test outcomes.

The 24-2 exam is the most common visual field test Carrot performs. In fact, Carrot users around the world conduct one of these exams every 12 seconds.



every 12 seconds a Carrot 24-2 exam is performed.

[24-2 Visual Field Exam Overview](#)

The 24-2 Visual Field Exam at a Glance

List of Ocular Diseases Monitored and Diagnoses Identified by the 24-2 Exam

Example 24-2 Reports

The 24-2C Visual Field Exam History

Billing and Coding for the 24-2 Eye Exam

Start Conducting the 24-2 Exam with Carrot

Clinical validation of Carrot's 24-2 exam



One study compared the Carrot head-mounted perimeter to the Humphrey Field Analyzer in glaucoma and optic nerve disease patients. Researchers described the output as “functionally indistinguishable” to the Humphrey Visual Field Analyzer.



Another study compared virtual reality visual field testing to Humphrey visual field testing in an academic ophthalmology practice. The results showed that virtual visual field testing had similar accuracy to Humphrey Carrot testing but offered additional advantages like shorter test duration, fewer false positives, and lower costs.



In 2024, Carrot received the Medical Device Single Audit Program (MDSAP) certification, verifying that our leading eye exam experience meets rigorous international standards, including ISO 13485:2016 and the stringent U.S. Food and Drug Administration (FDA) requirements.

30 days free.
No strings attached.

We are confident you'll love Carrot just like the 2,400+ doctors who have already made the switch.

[Start your 30-day trial](#)



24-2 Visual Field Exam Overview

[The 24-2 Visual Field Exam at a Glance](#)

List of Ocular Diseases Monitored and Diagnoses Identified by the 24-2 Exam

Example 24-2 Reports

The 24-2C Visual Field Exam History

Billing and Coding for the 24-2 Eye Exam

Start Conducting the 24-2 Exam with Carrot

The 24-2 Visual Field Exam at a Glance

Ophthalmologists appreciate the 24-2 visual field exam for its ability to detect subtle, early visual field defects and monitor progression in conditions like glaucoma and optic neuropathies, providing reliable, repeatable results crucial for managing chronic diseases.

However, they often express frustration with the test's limitations on traditional tabletop devices, like the lengthy and fatiguing process, reliance on patient cooperation, and the need for expert interpretation of complex reports. Thankfully, Carrot addresses some pain points by improving patient comfort and streamlining workflows.

Pros and cons of the 24-2 exam

We know the 24-2 exam is efficient, useful, and widely used, but there are pros and cons to incorporating this test into your routine.

PROS	CONS
This exam tests central and near-peripheral vision for a relatively comprehensive eye test.	Without proper patient preparation, the 24-2 exam can feel long and stressful.
You can identify subtle visual field defects, even before patients notice symptoms.	On traditional tabletop equipment, this exam can take several minutes per eye, which can impact clinic schedules and lead to patient fatigue.
Objective and repeatable results with standardized data support accurate disease monitoring.	Understanding reports and subtle changes in visual field requires expert interpretation.
As a cornerstone exam, the 24-2 is available on most visual field testing devices.	It's not a catch-all exam and more testing may be required to identify specific conditions.
The 24-2 is billable and reimbursable with patient insurance, including Medicare.	

24-2 Visual Field Exam Overview

The 24-2 Visual Field Exam at a Glance

[List of Ocular Diseases Monitored and Diagnoses Identified by the 24-2 Exam](#)

Example 24-2 Reports

The 24-2C Visual Field Exam History

Billing and Coding for the 24-2 Eye Exam

Start Conducting the 24-2 Exam with Carrot



11%

of Americans have diabetes

+25%

over a quarter of them will experience diabetic retinopathy

List of Ocular Diseases Monitored and Diagnoses Identified by the 24-2 Exam

Glaucoma	Glaucoma eye exams like the 24-2 monitor disease progression and assess the central visual field. Glaucoma can gradually affect peripheral vision, and the 24-2 test allows for early detection of functional changes in central and near-peripheral vision. It is often part of a comprehensive glaucoma workup for patients diagnosed with primary open-angle glaucoma, angle-closure glaucoma, or other related conditions.
Optic Neuropathy	Conditions like optic neuritis, optic neuropathy, or ischemic optic neuropathy can affect central vision. The 24-2 visual field exam is particularly useful in detecting and monitoring these conditions because it can identify visual field defects that may arise from optic nerve damage . This exam is crucial in tracking disease progression and assessing the impact of treatment in conditions that may compromise the optic nerve.
Neurological disorders	The 24-2 exam can help reveal neurological conditions that affect the visual processing centers in the brain. Conditions like stroke, brain tumors, multiple sclerosis, and traumatic brain injury may change or damage the brain's visual pathways , which the 24-2 test can help detect.
Macular Degeneration	Age-related macular degeneration and similar diseases usually affect the macula, but can also lead to peripheral visual field loss. The 24-2 test can detect these subtle changes, particularly in the early stages when functional changes may not be obvious on routine eye exams.
Diabetic Retinopathy	Diabetic retinopathy is a leading cause of blindness in adults. 11% of Americans have diabetes, and over a quarter of them will experience diabetic retinopathy. The 24-2 test can be used to detect early functional changes in the visual field, especially in the central region of vision that might not be detected by other diabetic eye exams.
Hypertension	Chronic hypertension can cause damage to the retinal blood vessels, and over time, this can lead to visual field changes. The 24-2 exam can help detect early signs of peripheral or central vision loss related to hypertension and identify potential retinal damage or other related conditions like hypertensive optic neuropathy.
Other Conditions	The 24-2 exam can also be used to monitor or diagnose: Retinal vein/artery occlusion • Retinal detachment • Retinitis pigmentosa • Autoimmune and inflammatory disorders • Tumors and brain lesions

24-2 Visual Field Exam Overview

The 24-2 Visual Field Exam at a Glance

List of Ocular Diseases Monitored and Diagnoses Identified by the 24-2 Exam

Example 24-2 Reports

The 24-2C Visual Field Exam History

Billing and Coding for the 24-2 Eye Exam

Start Conducting the 24-2 Exam with Carrot

Example 24-2 Reports

Carrot

RIGHT | Nov 7, 2024, 3:15 PM
Visual Field | Carrot

James Watson

7.4.2000 (24)

MRN
pa_ed3e842a

PATTERN
Central 24-2

STIMULUS
III, White

STRATEGY
BOLT

Rate of progression: $-0.7 \pm 5.4\%/yr$ (95% CI)
Not significant: Time between first and last exam must be >11 months

DATE	STRAT	DUR	FL	FP	FN	VFI	MD	PSD	GHT
4/6/24	BOLT	3:02	0%	0%	0%	99%	-5.49 db P < 5%	1.81 dB	Within Normal Limits
9/20/24	BOLT	2:59	0%	9%	0%	100%	-0.6 db	1.02 dB	Within Normal Limits
10/19/24	BOLT	2:56	0%	3%	0%	100%	-4.51 db	1.46 dB	Within Normal Limits
10/19/24	BOLT	3:45	0%	3%	0%	100%	-5.07 db P < 5%	1.52 dB	Within Normal Limits
10/20/24	BOLT	5:12	0%	0%	0%	98%	-5.71 db P < 5%	2.43 dB	Within Normal Limits
11/7/24	BOLT	3:11	0%	1%	0%	98%	-5.36 db P < 5%	2.03 dB	Within Normal Limits
11/7/24	BOLT	4:18	0%	4%	0%	98%	-5.36 db P < 5%	2.78 dB	Within Normal Limits

April 6, 2024

Grayscale

Threshold

Pattern Deviation

Baseline

24-2 Visual Field Exam Overview

The 24-2 Visual Field Exam at a Glance

List of Ocular Diseases Monitored and Diagnoses Identified by the 24-2 Exam


Example 24-2 Reports

The 24-2C Visual Field Exam History

Billing and Coding for the 24-2 Eye Exam

Start Conducting the 24-2 Exam with Carrot

Example 24-2 Reports



24-2 Visual Field Exam Overview

The 24-2 Visual Field Exam at a Glance

List of Ocular Diseases Monitored and Diagnoses Identified by the 24-2 Exam

Example 24-2 Reports

The 24-2C Visual Field Exam History

Billing and Coding for the 24-2 Eye Exam

Start Conducting the 24-2 Exam with Carrot

Example 24-2 Reports

LEFT | Nov 7, 2024, 3:15 PM
Visual Field | Carrot

James Watson

7.4.2000 (24)

MRN
pa_ed3e842a

STIMULUS
III, White

PATTERN
Central 24-2

STRATEGY
BOLT

DATE	STRAT	DUR	FL	FP	FN	VFI	MD	PSD	GHT
10/4/24	BOLT	3:37	0%	1%	0%	99%	-4.28 db	2.07 dB	Within Normal Limits
10/4/24	BOLT	2:43	0%	3%	0%	100%	-2.18 db	1.03 dB	Within Normal Limits
10/4/24	BOLT	2:57	0%	4%	0%	100%	-3.56 db	1.31 dB	Within Normal Limits
11/7/24	BOLT	3:19	0%	0%	0%	100%	-6.81 db P < 1%	2.06 dB	General Reduction of Sensitivity
11/7/24	BOLT	3:57	0%	1%	0%	100%	-4.68 db	1.78 dB	Within Normal Limits

24-2 Visual Field Exam Overview

The 24-2 Visual Field Exam at a Glance

List of Ocular Diseases Monitored and Diagnoses Identified by the 24-2 Exam


Example 24-2 Reports

The 24-2C Visual Field Exam History

Billing and Coding for the 24-2 Eye Exam

Start Conducting the 24-2 Exam with Carrot

Example 24-2 Reports



24-2 Visual Field Exam Overview

The 24-2 Visual Field Exam at a Glance

List of Ocular Diseases Monitored and Diagnoses Identified by the 24-2 Exam


Example 24-2 Reports

The 24-2C Visual Field Exam History

Billing and Coding for the 24-2 Eye Exam

Start Conducting the 24-2 Exam with Carrot

Example 24-2 Reports



24-2 Visual Field Exam Overview

The 24-2 Visual Field Exam at a Glance

List of Ocular Diseases Monitored and Diagnoses Identified by the 24-2 Exam

Example 24-2 Reports

The 24-2C Visual Field Exam History

Billing and Coding for the 24-2 Eye Exam

Start Conducting the 24-2 Exam with Carrot

Example 24-2 Reports

RIGHT | Apr 1, 2021, 8:00 AM
Visual Field | Carrot

John Doe

1.1.1951 (70)

MRN: pa_abcd1234
PATTERN: Central 24-2
STIMULUS: III, White
STRATEGY: Multiple

DATE	STRAT	DUR	FL	FP	FN	VFI	MD	PSD	GHT
4/1/19	BOLT	4:14	3/11	0%	5%	66%	-9.77 db P < 0.5%	7.87 dB P < 0.5%	Outside Normal Limits
2/1/20	BOLT	4:17	5/16	1%	6%	66%	-8.72 db P < 0.5%	9.19 dB P < 0.5%	Outside Normal Limits
4/1/21	BOLT	4:15	4/18	0%	12%	62%	-11.5 db P < 0.5%	8.37 dB P < 0.5%	Outside Normal Limits

April 1, 2019

Grayscale

Threshold

Threshold

Pattern Deviation

Pattern Deviation

Baseline

February 1, 2020

Grayscale

Threshold

Threshold

Pattern Deviation

Pattern Deviation

Baseline

24-2 Visual Field Exam Overview

The 24-2 Visual Field Exam at a Glance

List of Ocular Diseases Monitored and Diagnoses Identified by the 24-2 Exam

Example 24-2 Reports

The 24-2C Visual Field Exam History

Billing and Coding for the 24-2 Eye Exam

Start Conducting the 24-2 Exam with Carrot

Example 24-2 Reports



24-2 Visual Field Exam Overview

The 24-2 Visual Field Exam at a Glance

List of Ocular Diseases Monitored and Diagnoses Identified by the 24-2 Exam

Example 24-2 Reports

[The 24-2C Visual Field Exam History](#)

Billing and Coding for the 24-2 Eye Exam

Start Conducting the 24-2 Exam with Carrot

The 24-2C Visual Field Exam

The 24-2C visual field test was [developed in 2019](#) as a way to enhance the 24-2 exam and increase sensitivity in the central visual field.

Like its predecessor, the 24-2C divides the central 24 degrees of a patient’s vision into a grid to identify changes in the visual field. The 24-2C test includes the 54 points from the standard 24-2 test, plus 10 additional test points within the central 10 degrees of vision. The “C” refers to the increased sensitivity to detect central vision defects, which can point to glaucoma progression.

Before the introduction of this test, data with this level of detail required administering both the 24-2 and the 10-2 exam. [Combining the two](#) saves time, both for your workflow and your patients.

Academic references and clinical validation



[The introduction of the 24-2C visual field test](#) compares the 24-2 and 10-2 testing patterns, which traditionally were completed separately. Initial research shows that the 24-2C is comprehensive and can be treated as a hybrid test.



[Another study compared the performance](#) of the 24-2C and 10-2 test grids in patients with glaucoma and people suspected of having the condition. Both grids produced similar overall results, but the 10-2 detected more central visual field defects and showed better alignment with OCT scans. The 24-2C was faster, whereas the 10-2 provided more detailed insights.



[One study evaluated the effectiveness of 24-2C, 24-2, and 10-2 visual field tests](#) for detecting mild-stage glaucoma with central visual field defects. The 24-2C demonstrated superior detection accuracy, particularly in the upper-central visual field. It also significantly reduced test duration when compared to 24-2 and 10-2.

24-2 Visual Field Exam Overview

The 24-2 Visual Field Exam at a Glance

List of Ocular Diseases Monitored and Diagnoses Identified by the 24-2 Exam

Example 24-2 Reports

[The 24-2C Visual Field Exam History](#)

Billing and Coding for the 24-2 Eye Exam

Start Conducting the 24-2 Exam with Carrot

Conducting the 24-2C with Carrot

In 2024, Carrot added the 24-2C to our suite of visual field exams. This testing capability was automatically added to all Carrot devices, so there's no extra cost or need to upgrade. Simply update your software to the latest version to begin offering the 24-2C.

Whether you choose the traditional 24-2 or newer 24-2C, the testing and billing processes are the same.

30 days free.
No strings attached.

We are confident you'll love Carrot just like the 2,400+ doctors who have already made the switch.

[Start your 30-day trial](#)



24-2 Visual Field Exam Overview

The 24-2 Visual Field Exam at a Glance

List of Ocular Diseases Monitored and Diagnoses Identified by the 24-2 Exam

Example 24-2 Reports

The 24-2C Visual Field Exam History

[Billing and Coding for the 24-2 Eye Exam](#)

Start Conducting the 24-2 Exam with Carrot

Billing and Coding for the 24-2 Eye Exam

The 24-2 visual field exam is billable to insurance under CPT code 92083, which is designated for extended visual field exams.

According to the [Medicare Physician Fee Schedule \(MPFS\)](#), reimbursement usually ranges from \$40 to \$90 per test. Exact fees will, of course, depend on your location, setting, and payer-specific factors.

When is the 24-2 exam required

The 24-2 is a standard baseline exam for glaucoma patients. You can expect to re-examine people with early-stage glaucoma every six months to a year. For moderate to advanced glaucoma, cadence increases to every three to six months to monitor disease progression.

People often don't notice changes to their visual field, so this exam is also ordered for patients with diabetes, hypertension, or a family history of eye disease. People taking certain medications known to increase intraocular pressure may also need the 24-2 exam to monitor visual changes.

Is the 24-2 visual field test required for driver's licenses?

The 24-2 exam is more advanced than the basic field-of-vision checks required for driver's licensing. Currently, this exam isn't specifically mandated for drivers in the United States or Canada.

However, many states require a horizontal visual field of 120 degrees or better. In Canada, provinces have standards for visual acuity and peripheral vision, but don't specify the 24-2 test. That said, if you incorporate the 24-2 into your routine exams, the detailed results can support the information necessary for driver's licensing.

24-2 Visual Field Exam Overview

The 24-2 Visual Field Exam at a Glance

List of Ocular Diseases Monitored and Diagnoses Identified by the 24-2 Exam

Example 24-2 Reports

The 24-2C Visual Field Exam History

Billing and Coding for the 24-2 Eye Exam

[Start Conducting the 24-2 Exam with Carrot](#)

Start Conducting the 24-2 Exam with Carrot

The 24-2 visual field report is a wealth of diagnostic information when interpreted correctly. Diagnose, monitor, and manage a wide range of ocular and neurological conditions ranging from glaucoma to brain tumors. Carrot makes 24-2 testing easier, more comfortable, and less stressful — both for patients and technicians.

Download our complete Carrot Guide to Visual Exams to learn about all 13 tests you can perform with Carrot.

Ready to get started?

Schedule a demo or begin your 30-day free trial of Carrot to offer an exceptional eye exam experience in your practice.

[Book a demo](#)

[Start your 30-day trial](#)



Questions? Contact sales@carrot.io talk to a Carrot expert today.